

How to pair your hearing aids to your iPhone


Audiology Department

© Jan 2025
James Paget University Hospitals NHS
Foundation Trust
Review Date: Jan 2028
Code No version 1


Oticon Bluetooth hearing aids are 'Made For iPhone' (MFi) hearing aids, which means they can easily connect to your iPhone, iPad® or iPod touch®. You can then stream sound directly from your phone to your hearing aids, and use Oticon Companion to control your hearing aids from your screen. *Steps may vary depending on the iOS version installed on your device.




1. Go to Settings via your phone's home screen. Make sure Bluetooth is turned on




2. Scroll down a little then select Accessibility




3. From the options on the Accessibility menu, select Hearing Devices



4. At this point, restart your hearing aids and place them close to your device. They will remain in pairing mode for around 3 minutes



5. Once detected by your phone, tap the name of the Oticon hearing aid that you want to connect and confirm pairing in the pop-up. If you want to connect more than one hearing aid, repeat this step for the second hearing aid



6. Finally, when the devices have been connected, download Oticon Companion from the App Store®. Once installed and opened, it will automatically detect your paired hearing aids